# PeakTech®

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## PeakTech® 5310 A

**PH Messgerät /  
PH Meter**

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**PeakTech®**

## **1. Safety Precautions**

This product complies with the requirements of the following directives of the European Union for CE conformity: 2014/30/EU (electromagnetic compatibility), 2011/65/EU (RoHS). Pollution degree 2.

To ensure safe operation of the equipment, the following safety precautions must be observed.

Damages resulting from failure to observe these safety precautions are exempt from any legal claims whatever.

- \* Read these instructions carefully and make them accessible to subsequent users.
- \* Use this device only within its scope and specifications.
- \* Take all necessary safety precautions, especially when using acids.
- \* Wear personal safety equipment, gloves and safety glasses when handling chemical products.
- \* Never operate the device if it is not completely closed.
- \* Do not operate the device near strong magnetic fields (motors, transformers, etc.)
- \* Avoid strong vibrations of the device Before starting operation, the device should be stabilized to the ambient temperature (Important when transporting cold to warm rooms and vice versa)
- \* Do not make any technical changes to the device
- \* Only qualified service technicians may open the device, as well as perform maintenance and repair work

### **Cleaning the cabinet**

Clean only with a damp, soft cloth and a commercially available mild household cleanser. Ensure that no water gets inside the equipment to prevent possible shorts and damage to the equipment.

## 1.1. Features

- Digital bench PH meter, general purpose applications,
- easy operation and compact size
- Large 3 ½ - digit LCD display
- Calibration and adjustments at pH4 & pH7 with potentiometers
- Application: Education, School, Colleges, Laboratory, Industrial and Quality control

## 2. Specifications

Display	3 ½-digit, LCD-display max 1999 With backlight illumination
Measuring range	0.00 ~ 14.00 pH
Resolution	0,01 pH
Accuracy	+/- 0,1 pH
Measuring rate	1 – 2 / second
Calibration	External with potentiometer for pH4 and pH7
Operation Temperature	0°C...+50°C (32°F...122°C); < 85% RH
Automatic Temperature compensation	0°C...+50°C (32°F...122°C)
Power Supply	9V Battery (NEDA1604)
Dimensions (WxHxD)	128 x 68 x 25 mm
Weight	110g
Standard Accessories	Carrying Case, Manual, Battery, pH Probe, Screwdriver

### 3. Front Panel Description

1. PH-Electrode
2. containers designed for cleaning - or calibration buffer solution
3. BNC-socket for PH-Electrode
4. Battery compartment
5. Potentiometer for adjustment of PH-7
6. Potentiometer for adjustment of PH-4
7. Power switch
8. Display

## **4. Measurement method**

### **4.1. pH measurement**

After calibrating the measuring device and the pH electrode, the device is ready for measurement.

- 1.) Connect the electrode to the BNC connector of the device.
- 2.) Remove the protective cap from the electrode.
- 3.) Switch the measuring device on by pressing the ON / OFF button.
- 4.) Immerse the pH electrode in the solution to be measured.
- 5.) Slowly stir the electrode in the solution until you get a stable reading.
- 6.) After you have taken a measurement, rinse the electrode with distilled water!
- 7.) Put the protective cap back on the electrode.

## **4.2. Calibration**

For calibration, a standard buffer solution pH6.86, pH4.01 and pH 9.18 is required, which can be obtained from specialist dealers.

- 1.) Immerse the electrode in a standard phosphate buffer solution with a pH value of 6.86 (at a temperature of 25 ° C) and stir gently.
- 2.) Turn the trimmer pH7 with a screwdriver until the measured value corresponds to 6.86pH.
- 3.) Clean the electrode.
- 4.) Immerse the electrode in a standard buffer solution of pH 4.01 (at a temperature of 25 ° C). Turn the pH4 trimmer with a screwdriver until the measured value equals 4.01pH.
- 5.) Clean the electrode.
- 6.) Control: Immerse the electrode in a standard buffer solution of pH 9.18 (below a temperature of 25 ° C) and stir gently until the display value is stable. If the measured value is within the acceptable error range, the calibration is complete. If not, repeat point 4.

### **Important:**

The pH range of the instrument needs to be recalibrated if:

- The electrode has not been used for a long time since the last calibration.
- The device has been used frequently and extensively.
- When the highest possible accuracy is required.
- After the electrode has been replaced with a new electrode.

## **6. Replacement of Battery**

- 1) Remove the battery cover on the rear side of the device.
- 2) Take out the old battery and insert new 9V-battery (9V NEDA 1604 or equivalent) according to the polarity.
- 3) Replace the battery compartment cover back on.

## **Statutory Notification about the Battery Regulations**

The delivery of many devices includes batteries, which for example serve to operate the remote control. There also could be batteries or accumulators built into the device itself. In connection with the sale of these batteries or accumulators, we are obliged under the Battery Regulations to notify our customers of the following:

Please dispose of old batteries at a council collection point or return them to a local shop at no cost. The disposal in domestic refuse is strictly forbidden according to the Battery Regulations. You can return used batteries obtained from us at no charge at the address on the last side in this manual or by posting with sufficient stamps.

Batteries, which contain harmful substances, are marked with the symbol of a crossed-out waste bin, similar to the illustration shown left. Under the waste bin symbol is the chemical symbol for the harmful substance, e.g. „Cd“ for cadmium, „Pb“ stands for lead and „Hg“ for mercury.

You can obtain further information about the Battery Regulations from the [Bundesministerium für Umwelt, Naturschutz und Reaktorsicherheit](#) (*Federal Ministry of Environment, Nature Conservation and Reactor Safety*).

*All rights, also for translation, reprinting and copy of this manual or parts are reserved. Reproductions of all kinds (photocopy, microfilm or other) only by written permission of the publisher.*

*This manual is according the latest technical knowing. Technical alterations reserved.*

*Misprints and errors are reserved.*

*We herewith confirm that the units are calibrated by the factory according to the specifications as per the technical specifications.*

*We recommend to calibrate the unit again, after 1 year.*

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© **PeakTech**® 12/2021 Ehr.
PeakTech Prüf- und Messtechnik GmbH